



## OP-SYSTEM

### Veterinary Surgery and Treatment Table

- High-end veterinary special-surgery table for ophthalmology and advanced orthopedics
- Developed for complex procedures and surgery
- Frame of maximum stability
- Powder-coated according to RAL 7035
- Surgery also possible while seating
- Electrical height adjustment from 59 - 99 cm
- Optionally moveable by double castors or fixed table
- Height adjustment by hand switch
- Lateral up to 40° electric tilting
- Longitudinal up to 15° electric tilting in both directions
- X-ray translucent plate
- Optionally special version with sliding top for ophthalmology procedures or C-bow x-ray application
- Castors with directional locking and central locking mechanism by actuating a single pedal

#### Basic equipment:

- OP-System with castors
- OP rails (2 pieces lengthwise with fixed top; 4 pieces lengthwise with sliding top)
- Quick release tie-down cleats (4 pieces)
- Hand switch
- Optional:* Special version with sliding top
- Optional:* Joystick



Special version with sliding top

## OP-SYSTEM

moveable by double castors with blocking mechanism



### Electrical height adjustable (power cord)

Top	Item Code
fixed top	652. 521
sliding top	652. 541

### Battery-powered height adjustable

Top	Item Code
fixed top	672. 521
sliding top	672. 541

## OP-SYSTEM

fixed

### Electrical height adjustable (power cord)

Top	Item Code
fixed top	652. 520
sliding top	652. 540

### Battery-powered height adjustable

Top	Item Code
fixed top	672. 520
sliding top	672. 540

## Accessoires for OP-SYSTEM


### Available accessory items

Item	Item Code
Quick-release tie-down cleat*	602. 601
Infusion pole*	602. 602
Adjustable side supports (pair)*	602. 604
Narcotic bow*	602. 605
Leg support with tie-down cleat*	602. 606
Tilting stainless steel dish with grid*	602. 607
Rail extension (200 mm, pair)*	602. 610
Instrument tray (300x200 mm)*	602. 611
Arm rest for ophthalmology*	602. 614
Top extension (500x300 mm)**	602. 615
Top extension (300x300 mm)**	602. 616
Insert for OP-System**,***	602. 635

\* suitable for OP rails

\*\* only compatible with OP-System

\*\*\* only available for sliding top version



Top extension (500x300 mm)  
**602.601**

Top extension (500x300 mm)  
**602.601**



Insert for OP-System  
**602.635**



Power supply: 230V 50/60Hz 3A

Lifting capacity: 120 kg