




VET LIFT

Veterinary Surgery and Treatment Table

- High-end scissor type surgery and treatment table with versatile use in veterinary medicine
- Frame of maximum stability; powder coated according to RAL 7035
- Electrical, battery-powered or hydraulic height adjustment by foot switch from 31-101 cm
- Plate (stainless steel or synthetic top) with edge and drain (130 x 60 cm) as well as 14 tie-down hooks
- Optionally fixed with 2 castors installed on one side (allow smooth movement of the table by lifting on the opposite side) or moveable by 4 double castors with blocking mechanism
- Top optionally fixed or tiltable (continuously adjustable tilt angle of the top (up to 11°) with comfortable one-hand operation supported by gas spring)



Special version with integrated scale

- 0-100 kg in 50 g steps
- Animal weighing function (averaging method)
- Battery-operated (220 hours working hours)
- Automatic shutdown after 3 minutes
- Delivery without batteries (6 pieces AA)



Special version with integrated scale

VET LIFT

fixed with small castors on one side



Hydraulic height adjustable

Top	Item Code
fixed stainless steel top	769. 605
tilting stainless steel top	769. 606
fixed synthetic top	769. 615
tilting synthetic top	769. 616

Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759. 605
tilting stainless steel top	759. 606
fixed synthetic top	759. 615
tilting synthetic top	759. 616

Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779. 605
tilting stainless steel top	779. 606
fixed synthetic top	779. 615
tilting synthetic top	779. 616

VET LIFT

movable by 4 double castors (with blocking mechanism)



Hydraulic height adjustable

Top	Item Code
fixed stainless steel top	769. 601
tilting stainless steel top	769. 602
fixed synthetic top	769. 611
tilting synthetic top	769. 612

Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759. 601
tilting stainless steel top	759. 602
fixed synthetic top	759. 611
tilting synthetic top	759. 612

Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779. 601
tilting stainless steel top	779. 602
fixed synthetic top	779. 611
tilting synthetic top	779. 612

VET LIFT with scale

fixed with small castors on one side



Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	759. 607
fixed synthetic top	759. 617

Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	779. 607
fixed synthetic top	779. 617

