



ECO LIFT

Veterinary Surgery and Treatment Table

- Chassis made of high-end steel and aluminium
- Powder-coated surface according to RAL card
- Optionally electric (power cord) or hydraulic height adjustment from 72 - 102 cm
- Frame of maximum stability
- Optionally stainless steel or synthetical table top plate with drain (130 x 60 cm)
- Solid stand of the table due to stable and high-quality base frame
- Continuously adjustable tilting angle of the table top up to almost vertical position (approx. 90°)
- Special shape of the base frame enables power supply via a floor socket or power column
- Chassis available in five different colours



Special version with heated stainless steel top



Special version with washing table

ECO LIFT

fixed



Electrical height adjustable (power cord)

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top	651. 600	. xxxx
tilting synthetic top	651. 610	. xxxx

Hydraulic height adjustable

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top	661. 600	. xxxx
tilting synthetic top	661. 610	. xxxx

ECO LIFT Heating

fixed



Electrical height adjustable (power cord)

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top <i>heated!</i>	651. 600	. xxxx . SH

Hydraulic height adjustable

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top <i>heated!</i>	661. 600	. xxxx . SH

ECO LIFT Washing Table

fixed



Electrical height adjustable (power cord)

Top	Item Code	Colour Code (Chassis)
Washing Table Stainless Steel	651. 600	. xxxx - W

Hydraulic height adjustable

Top	Item Code	Colour Code (Chassis)
Washing Table Stainless Steel	661. 600	. xxxx - W

Colour card



red-violet (Nr. 4002)



signal blue (Nr. 5005)



daffodil yellow (Nr. 1007)



light grey (Nr. 7035)



light green (Nr. 6027)



Power supply: 230V 50Hz 0,63A

Lifting capacity: 120 kg