




AERON

Surgery table for C-arm application

- Electrical height adjustable from 80 - 110 cm
- Longitudinal 15° electrical tilting
- Lateral 20° electrical tilting
- X-ray translucent plate (exact area dependant on chosen version)
- Suitable for arthroscopy in small animals (extensive accessories available on request)
- Battery-operated or with power cord
- Movable on double castors
- Fixation on lockable, anti-slip rubber feet
- Lifting capacity: 120 kg

Basic equipment

- OP rails
- 4 pieces quick-release tie-down cleats
- Joystick remote control



Lateral 20° electric tilt

Version with V-Top

- Centrally divided frame
- One frame half can be removed at a time

AERON

movable on castors with blocking mechanism



V-Top

Type	Item Code
Battery-powered	776.002
Electric with power cord	756.002



Flat Top

Type	Item Code
Battery-powered	776.522
Electric with power cord	756.522

Attention: Frame cannot be removed with flat top version of the table!
Basic equipment without shown side supports
(available as accessory)

Accessoires for AERON

Item	Item Code
Quick-release tie-down cleats	602.601
Infusion pole	602.602
Adjustable side supports (pair)	602.604
Narcotic bow	602.605
Leg support with tie-down cleat	602.606
Tilting stainless steel dish with grid	602.607
Rail extension (200 mm, pair)	602.610
Instrument tray (300 x 200 mm)	602.611
Arm rests for ophtalmology	602.614
Top extension 500 x 300 mm	602.615
Top extension 300 x 300 mm	602.616



Power supply: 230V 50/60Hz max. 3A

Loading capacity: 120 kg