



ECO VET electrical Veterinary Surgery and Treatment Table

- Electrical height adjustment from 54 - 100 cm
- Optionally high-quality stainless steel or synthetic table top with edge and drain (130 x 60 cm)
- 14 tie-down cleats
- Optionally battery-operated or electrical by power cord height adjustment
- Optionally stainless steel or synthetic table top plate; available in fixed or tilting edition
- Foot switch for height adjustment



Seamless plate
Completely lined underneath

ECO VET

fixed



Electrical height adjustable (power cord)

Top	Item Code
fixed stainless steel top	758.605
tilting stainless steel top	758.600
fixed synthetic top	758.615
tilting synthetic top	758.610

Battery-powered height adjustable

Top	Item Code
fixed stainless steel top	778.605
tilting stainless steel top	778.600
fixed synthetic top	778.615
tilting synthetic top	778.610

Accessories for ECO VET

Accessory:

Item	Item Code
OP rails (set of 4 pieces)	602.600
Quick-release tie-down cleat*	602.601
Infusion pole*	602.602
Adjustable side support (pair)*	602.604
Narcotic bow*	602.605
Leg support with tie-down cleat*	602.606
Tilting stainless steel dish with grid*	602.607
Rail extension (200 mm / pair)*	602.610
Instrument tray (300x200 mm)*	602.611
Arm rest for ophthalmology*	602.614

* suitable for OP rails



Power supply: 230V 50Hz 1,5A

Lifting capacity: 120 kg