



TREND

Veterinary Surgery and Treatment Table

- Chassis made of high-end steel and aluminium
- Powder-coated surface according to RAL card
- Electric (power cord) height adjustment from 72 - 102 cm
- Operation via foot switch on both sides
- Frame of maximum stability
- Optionally stainless steel or synthetical table top plate with drain (130 x 60 cm)
- Equipped with castors on one side (allow smooth movement of the table by lifting on the opposite side)
- Continuously adjustable tilting angle of the table top up to almost vertical position (approx. 90°)
- Comfortable one-hand operation supported by gas spring
- Special shape of the base frame enables power supply via a floor socket or power column
- Available in five different colours



Special version with heated stainless steel top

TREND

moveable with castors on one side



Electrical height adjustable (power cord)

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top	754.601	.xxxx
tilting synthetic top	754.611	.xxxx

TREND Heating

moveable with castors on one side



Electrical height adjustable (power cord)

Top	Item Code	Colour Card (Chassis)
tilting stainless steel top <i>heated!</i>	754.601	.xxxx .SH

Colour card



red-violet (Nr. 4002)



daffodil yellow (Nr. 1007)



light green (Nr. 6027)



signal blue (Nr. 5005)



light grey (Nr. 7035)



Power supply: 230V 50Hz 0,63A

Lifting capacity: 120 kg