



ECO LIFT V-TOP

Veterinary Surgery and Treatment Table

- Chassis made of high-end steel and aluminium
- Powder-coated surface according to RAL card
- Optionally electric (power cord) or hydraulic height adjustment from 72 - 102 cm
- Frame of maximum stability
- Removeable, high-end stainless steel table top (total size of 143x50 cm)
- Solid standing of the table due to robust base frame
- Mechanically adjustable tilt angle of the table top (up to 45°)
- Mechanically variable V-TOP (up to 90°) for comfortable adjustment
- Including drain gutter and drip tray (extendable)
- Including 2 rail (on the length side)
- Chassis available in five different colours



ECO LIFT V-TOP

fixed



Electrical height adjustable (power cord)

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top	651. 531	. xxxx

Hydraulic height adjustable

Top	Item Code	Colour Code (Chassis)
tilting stainless steel top	661. 531	. xxxx

Colour card



red-violet (Nr. **4002**)



daffodil yellow (Nr. **1007**)



light green (Nr. **6027**)



signal blue (Nr. **5005**)



light grey (Nr. **7035**)



Power supply: 230V 50Hz 0,63A

Lifting capacity: 120 kg