



OP-System V-TOP

Veterinary Surgery and Treatment Table

- High-end veterinary special-surgery table for ophthalmology and advanced orthopedics
- Developed for complex procedures and surgery
- Surgery also possible while seating
- Electrical height adjustment from 58 - 98 cm
- Moveable by double castors with blocking mechanism
- Lateral up to 40° electric tilting
- Continuously variable V-Top position (possible up to 90°) with comfortable one-hand operation supported by gas spring
- Longitudinal up to 10° electric tilting in both directions
- Castors with directional locking and central locking mechanism by actuating a single pedal

Basic equipment:

OP-System with castors

OP rails (2 pieces)

Hand switch

Optional: Joystick



OP-SYSTEM V-TOP

Standard



Electrical height adjustable (power cord)

Movement	Item Code
not moveable (fixed)	652. 530
moveable by double castors with blocking mechanism	652. 531

Battery-powered height adjustable

Movement	Item Code
not moveable (fixed)	672. 530
moveable by double castors with blocking mechanism	672. 531



Power supply: 230V 50/60Hz 3A

Lifting capacity: 120 kg